

Product Datasheet

ZO1 antibody (orb556544)

Description

Rabbit polyclonal antibody to ZO1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the C-terminus region of human ZO-1. The exact sequence is proprietary.

Target

ZO1

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HCT116 , 293T , A431 , HeLa

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


195

Uniprot ID
[Q07157](https://www.uniprot.org/entry/Q07157)
Dilution Range


WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of 293T whole cell...



Western blot analysis of HCT116 cell...