

Product Datasheet

RAP2C, member of RAS oncogene family Antibody (orb556546)

Description

Rabbit polyclonal antibody to RAP2C



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunohistochemical staining of N87 xeno...

Immunogen

Full length human RAP2C Recombinant protein.

Target


RAP2C

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.




Western blot analysis of MOLT4 cells usi...

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunohistochemical staining of mouse du...

Note

For research use only

Application notes

Positive control: Molt-4

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

21

Uniprot ID

[Q9Y3L5](https://www.uniprot.org/entry/Q9Y3L5)

Dilution Range

WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.