


Product Datasheet

SLC25A13 antibody (orb556579)

Description

Rabbit polyclonal antibody to SLC25A13


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of HeLa
cell...

Immunogen

Recombinant protein encompassing a sequence within the center region of human SLC25A13. The exact sequence is proprietary.

Target


SLC25A13

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.



Immunohistochemical
staining of N87
xeno...

Concentration

0.61 mg/ml

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse liver , rat liver , 293T , A431 , HeLa , HepG2



Immunohistochemical
staining of mouse
ki...

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

74

Uniprot ID
[Q9UJS0](https://www.uniprot.org/entry/Q9UJS0)
Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.