

Product Datasheet

**aldo-keto reductase family 1 member B10
Antibody (orb556620)**

Description

Rabbit polyclonal antibody to AKR1B10



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation


Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of HeLa
cell...



Immunohistochemical
staining of OVCA3
X...

Immunogen

Full length human AKR1B10 Recombinant protein.

Target

AKR1B10

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: mouse colon , rat colon

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

36

Uniprot ID


[O60218](https://www.uniprot.org/entry/O60218)

Dilution Range

WB: 1:500-1:10000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis
of MOLT4(Lane 1), ...