


## Product Datasheet

### PCAF antibody (orb556629)


## Description

Rabbit polyclonal antibody to PCAF

|                            |                                                                                                                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                               |
| <b>Tested Applications</b> | WB                                                                                                                                                                                                                         |
| <b>Immunogen</b>           | Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human PCAF. The exact sequence is proprietary.                                                                        |
| <b>Target</b>              | PCAF                                                                                                                                                                                                                       |
| <b>Preservatives</b>       | 0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol                                                                                                                                                                                |
| <b>Form/Appearance</b>     | Liquid: 1XPBS, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.                                                                                                                                           |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                                                                                                    |
| <b>Storage</b>             | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                      |
| <b>Application notes</b>   | Positive control: Mouse brain , 293T , A431 , A375 , HeLa nuclear extract                                                                                                                                                  |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                 |
| <b>Purity</b>              | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| <b>MW</b>                  | 93                                                                                                                                                                                                                         |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q92831">Q92831</a>                                                                                                                                                                  |
| <b>Dilution Range</b>      | WB: 1:500-1:3000                                                                                                                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                            |



Western blot analysis of 293T whole cell...



Western blot analysis of mouse brain lysate...