

Product Datasheet

AUTS2 antibody (orb556638)

Description

Rabbit polyclonal antibody to AUTS2

Species/Host

Rabbit

Reactivity

Human, Mouse, Zebrafish

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human AUTS2. The exact sequence is proprietary.

Target

AUTS2

Preservatives

0.025% ProClin 300, 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: NIH-3T3, zebrafish brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


139

Uniprot ID
[Q8WXX7](#)
Dilution Range


WB: 1:500-1:3000

Expiration Date


12 months from date of receipt.



Western blot analysis of Zebrafish tissue...



Western blot analysis of 293T cells using...



Western blot analysis of NIH-3T3 cells using...