


## Product Datasheet

### Fatty Acid Synthase antibody (orb556642)

## Description

Rabbit polyclonal antibody to Fatty Acid



## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat


## Conjugation

Unconjugated

## Tested

## Applications

ICC, IHC-Fr, IHC-P, IP, WB



Immunofluorescence  
analysis of HeLa  
cell...

## Immunogen


Recombinant protein encompassing a sequence within the N-terminus region of human Fatty Acid Synthase. The exact sequence is proprietary.

## Target

Fatty Acid Synthase

## Preservatives

0.025% ProClin 300, 1XPBS (pH7), 1% rAlbumin, 20% Glycerol



Immunohistochemical  
staining of human  
co...

## Form/Appearance


Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

## Concentration

1 mg/ml

## Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunohistochemical  
staining of human  
br...

## Note

For research use only

## Application notes

Positive control: Mouse brain, Molt-4, PC-12, Rat2, 293T

## Isotype

IgG

## Clonality

Polyclonal

## Purity

Purified by antigen-affinity chromatography.

## MW

273

## Uniprot ID

[P49327](https://www.uniprot.org/entry/P49327)

## Dilution Range

WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000, IP: 1:100-1:1000

## Expiration Date

12 months from date of receipt.