

Product Datasheet

58K Golgi protein antibody (orb556656)

Description

Rabbit polyclonal antibody to 58K Golgi protein



Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human 58K Golgi protein. The exact sequence is proprietary.

Target

58K Golgi protein

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.98 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse liver , HepG2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

59

Uniprot ID

[O95954](#)


Dilution Range

WB: 1:500-1:3000 , IHC-P: 1:100-1:1000


Expiration Date

12 months from date of receipt.

Immunohistochemical staining of human he...



Western blot analysis of Hep G2 cells us...



Western blot analysis of mouse liver usi...