

Product Datasheet

MCM7 antibody (orb556685)

Description

Rabbit polyclonal antibody to MCM7

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ChIP, ICC, IHC-P, IP, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human MCM7. The exact sequence is proprietary.

Target

MCM7

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.76 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: H1299

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

81

Uniprot ID


[P33993](#)

Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000, IP: 1:100-1:500

Expiration Date


12 months from date of receipt.



Chromatography analysis of HeLa cell
usi...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of Cal27 Xe...