


Product Datasheet

Fatty Acid Synthase antibody (orb556692)

Description

Rabbit polyclonal antibody to Fatty Acid



Species/Host

Rabbit

Reactivity

Human, Mouse


Conjugation

Unconjugated

Tested

Applications

ICC, IHC-P, IP, WB



Immunofluorescence
analysis of HeLa
cell...

Immunogen

Recombinant protein encompassing a sequence within the N-terminus region of human Fatty Acid Synthase. The exact sequence is proprietary.

Target


Fatty Acid Synthase

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.




Immunofluorescence
analysis of HeLa
cell...

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunohistochemical
staining of OECM1
xe...

Note

For research use only

Application notes

Positive control: Molt-4 , 293T , A431 , HeLa , HepG2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

273

Uniprot ID

[P49327](#)

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IP: 1:100-1:500

Expiration Date

12 months from date of receipt.