

Product Datasheet

NFKB inhibitor alpha Antibody (orb556701)

Description

Rabbit polyclonal antibody to IKB alpha

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine

Conjugation

Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human IKB alpha. The exact sequence is proprietary.

Target

IKB alpha

Preservatives

0.025% ProClin 300, 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.

Concentration

0.16 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: BCL-1, C2C12, NIH-3T3, HT-29, 293T, A431, HeLa, HepG2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

36

Uniprot ID


[P25963](#)

Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of mouse ki...



Western blot analysis of Various whole c...