

Product Datasheet

NUF2 component of NDC80 kinetochore complex
Antibody (orb556711)

Description

Rabbit polyclonal antibody to NUF2


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ICC, WB

Applications

Immunofluorescence
analysis of A431
cell...

Immunogen

Recombinant protein corresponding to human CDCA1

Target


NUF2

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.



Immunofluorescence
analysis of U2OS
cell...

Concentration

0.35 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


54

Uniprot ID
[Q9BZD4](https://www.uniprot.org/entry/Q9BZD4)
Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Western blot
analysis of HeLa
cells usin...