

Product Datasheet

F-box and leucine rich repeat protein 3 Antibody (orb556714)

Description

Rabbit polyclonal antibody to FBXL3



Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human FBXL3. The exact sequence is proprietary.

Target

FBXL3

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.41 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: C2C12 , Molt-4

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

49

Uniprot ID


[Q9UKT7](https://www.uniprot.org/entry/Q9UKT7)

Dilution Range


WB: 1:500-1:3000 , IF: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of A431 cell...



Western blot analysis of Molt-4 cells us...



Western blot analysis of C2C12 whole cel...