


Product Datasheet

CYFIP2 antibody (orb556719)

Description	Rabbit polyclonal antibody to CYFIP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human CYFIP2. The exact sequence is proprietary.
Target	CYFIP2
Preservatives	0.025% ProClin 300. 1XPBS pH 7, 20% Glycerol
Form/Appearance	Liquid: 1XPBS, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	1 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	148
Uniprot ID	Q96F07
Dilution Range	WB: 1:5000-1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of Various whole c...