


Product Datasheet

CaMKII delta antibody (orb556746)

Description

Rabbit polyclonal antibody to CaMKII delta


Species/Host

Rabbit

Reactivity

Human, Mouse, Rabbit, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications

Immunofluorescence
analysis of HCT116
ce...

Immunogen

Recombinant protein encompassing a sequence within the center region of human CaMKII delta. The exact sequence is proprietary.

Target


CaMKII delta

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.


Concentration

0.41 mg/ml

Immunohistochemical
staining of rat hear...

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: A549 , H1299 , HCT116 , 293T , A431 , HeLa , HepG2 , mouse brain , rat brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


56

Uniprot ID
Q13557
Dilution Range

WB: 1:500-1:50000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IP: 1:100-1:500

Expiration Date

12 months from date of receipt.


Immunohistochemical
staining of mouse
mu...