


Product Datasheet

Fructose 6 Phosphate Kinase antibody (orb556749)


Description

Rabbit polyclonal antibody to Fructose 6 Phosphate

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human Fructose 6 Phosphate Kinase. The exact sequence is proprietary.
Target	Fructose 6 Phosphate Kinase
Preservatives	0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	0.92 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: mouse brain , rat brain , HeLa
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	85
Uniprot ID	P08237
Dilution Range	WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell extr...



Western blot analysis of various tissue ...