

Product Datasheet

phosphodiesterase 1A Antibody (orb556763)

Description

Rabbit polyclonal antibody to PDE1A

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human PDE1A. The exact sequence is proprietary.

Target

PDE1A

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.74 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: H1299, HeLa, mouse brain, rat brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

61

Uniprot ID


[P54750](#)

Dilution Range


WB: 1:500-1:3000, IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of MDAMB468...



Western blot analysis of H1299(Lane 1), ...



Western blot analysis of H1299(Lane 1), ...