

Product Datasheet

enolase 3 Antibody (orb556851)

Description

Rabbit polyclonal antibody to ENO3

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human ENO3. The exact sequence is proprietary.

Target

ENO3

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.46 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HepG2 , Molt-4 , Raji , Mouse brain , Rat2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

47

Uniprot ID


[P13929](#)

Dilution Range


WB: 1:1000-1:10000 , IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human ga...



Western blot analysis of Hep G2(Lane 1),...



Western blot analysis of Hep G2(Lane 1),...