


Product Datasheet

heterogeneous nuclear ribonucleoprotein C (C1/C2) Antibody (orb556856)

Description

Rabbit polyclonal antibody to hnRNP



Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IHC-P, IP, WB



Immunofluorescence analysis of HeLa cell...

Immunogen


Recombinant protein encompassing a sequence within the center region of human hnRNP C1/C2. The exact sequence is proprietary.

Target

hnRNP C1/C2

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol



Immunohistochemical staining of mouse br...

Form/Appearance


Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.79 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunohistochemical staining of Cal27 Xe...

Note

For research use only

Application notes

Positive control: Raji , Neuro2A , NIH-3T3 , BCL-1 , Raw264.7 , C2C12 , PC-12 , Rat2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

34

Uniprot ID

P07910

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.