

Product Datasheet

SSA1 antibody (orb556870)

Description

Rabbit polyclonal antibody to SSA1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human SSA1. The exact sequence is proprietary.

Target

SSA1

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.

Concentration

1.83 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: TRIM21-transfected 293T

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

54

Uniprot ID

[P19474](#)

Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000, IP: 1:100-1:500

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of A431 cell...



Immunohistochemical staining of DLD-1 cell...



Immunoprecipitation analysis of 293T cell...