

Product Datasheet

**tripartite motif containing 32 Antibody
(orb556889)**

Description

Rabbit polyclonal antibody to TRIM32


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human TRIM32. The exact sequence is proprietary.

Target

TRIM32

Preservatives

0.01% Thimerosal, 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: DDDDK-human TRIM32-transfected 293T cells

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


72

Uniprot ID
[Q13049](https://www.uniprot.org/entry/Q13049)
Dilution Range


WB: 1:500-1:20000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of A375 cell...



Immunohistochemical staining of mouse br...



Western blot analysis of 293T whole cell...