

Product Datasheet

Activin Receptor Type IIA antibody (orb556898)

Description

Rabbit polyclonal antibody to Activin Receptor Type

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human Activin Receptor type IIA. The exact sequence is proprietary.

Target

Activin Receptor Type IIA

Preservatives

0.025% ProClin 300. 1XPBS pH 7, 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.86 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

58

Uniprot ID

P27037

Dilution Range

WB: 1:5000-1:20000

Expiration Date

12 months from date of receipt.

