



Product Datasheet

UGT1A antibody (orb556901)


Description	Rabbit polyclonal antibody to UGT1A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human UGT1A. The exact sequence is proprietary.
Target	UGT1A
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 20% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	0.47 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: mouse liver , 293T , A431 , HepG2
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
Dilution Range	WB: 1:500-1:10000 , IHC-P: 1:100-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human he...



Western blot analysis of Various whole c...



Western blot analysis of Mouse tissue us...