



Product Datasheet

GNAT2 antibody (orb556912)

Description	Rabbit polyclonal antibody to GNAT2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Recombinant protein encompassing a sequence within the center region of human GNAT2. The exact sequence is proprietary.
Target	GNAT2
Preservatives	0.025% ProClin 300. 1XPBS pH 7, 20% Glycerol
Form/Appearance	Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	1 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	40
Uniprot ID	P19087
Dilution Range	IHC-Fr: 1:100-1:1000
Expiration Date	12 months from date of receipt.



Immunocytochemistry
analysis of mouse
re...



Immunocytochemistry
analysis of mouse
re...