


Product Datasheet

Histone H3 antibody (orb556958)

Description

Rabbit polyclonal antibody to Histone H3


Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of HeLa
cell...

Immunogen

Full length human Histone H3.3B Recombinant protein.

Target


Histone H3

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.


Immunohistochemical
staining of mouse
sp...

Concentration

0.16 mg/ml

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse brain , rat brain , 293T , A431 , HeLa , HepG2


Immunohistochemical
staining of mouse
lu...

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

15

Uniprot ID

P84243

Dilution Range

WB: 1:500-1:20000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.