

Product Datasheet

microtubule associated protein tau Antibody (orb556964)

Description

Rabbit polyclonal antibody to Tau



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ICC, WB

Applications

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the N-terminus region of human Tau. The exact sequence is proprietary.

Target

Tau

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1.03 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: SK-N-SH , IMR32 , SK-N-AS

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

79

Uniprot ID


[P10636](#)

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000


Expiration Date

12 months from date of receipt.



Immunofluorescence analysis of U-87 MG

C...



Western blot analysis of Various whole c...