

Product Datasheet

Mucin 2 antibody (orb557085)

Description	Mouse monoclonal antibody to Mucin 2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	Contains 0.09% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,IHC-P (1:10-1:40)
Expiration Date	12 months from date of receipt.