

Product Datasheet

Tubulin (Alpha) antibody (orb557115)

Description	Mouse monoclonal antibody to Tubulin (Alpha)
Species/Host	Mouse
Reactivity	Human, Mouse, Other, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgM-K
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.