

Product Datasheet

EBV-LMP1 antibody (orb557161)

Description	Mouse monoclonal antibody to EBV-LMP1
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, IP, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.