


Product Datasheet

HPV-16 E6 antibody (orb557165)

Description	Mouse monoclonal antibody to HPV-16 E6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a-K
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HPV-16 E6...



Immunohistochemical staining of HPV-16 E6...