

Product Datasheet

CKS1B antibody (orb55717)

Description	Rabbit polyclonal antibody to CKS1B.
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Goat, Guinea pig, Human, Mouse, Rat, Yeast, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	Cks1b
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EFEYRHVMLPKDIAKLVPKTHLMSESEWRNLGVQSQGWVHYMIHEPEPH
Clonality	Polyclonal
MW	10kDa
Uniprot ID	P61025
NCBI	NP_058600
Expiration Date	12 months from date of receipt.