

Product Datasheet

CD54 (ICAM-1) antibody (orb557221)

Description	Mouse monoclonal antibody to CD54 (ICAM-1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a-K
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.