

Product Datasheet

CD66e (CEACAM5) antibody (orb557229)

Description	Mouse monoclonal antibody to CD66e (CEACAM5)
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1-K
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.