

Product Datasheet

MARC2 antibody (orb55725)

Description	Rabbit polyclonal antibody to MARC2 .
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human MARC2
Target	MARC2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ACPRCILTTPDTPDTGVDRKQPLDTLKSRYRLCDPSERELYKLSPLFGIYY
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q969Z3
NCBI	NP_060368
Expiration Date	12 months from date of receipt.