

Product Datasheet

ANTXR1 (phospho-Tyr382) antibody (orb5578)

Description	ANTXR1 (phospho-Tyr382) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ANTXR1 around the phosphorylation site of Tyr382 (AS(p-Y)YG)
Target	phospho-TEM8 (Tyr382)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	Q9H6X2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.