

Product Datasheet

Human beta-1 protein (orb5581)

Description	Recombinant of human beta-1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 1Mnacl,0.05mMHeparin,50mMsodium acetate ,containing 0.01% sodium azide, pH 5.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Protein Sequence	Ile31-Pro275
Source	Human Cells
Uniprot ID	Q86TM8
Expiration Date	6 months from date of receipt.