



Product Datasheet

Tuberin (phospho-Thr1462) antibody (orb5583)

Description	Tuberin (phospho-Thr1462) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tuberin around the phosphorylation site of Thr1462 (GY(p-T)IS)
Target	Phospho-Tuberin (Thr1462)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	200 kDa
Uniprot ID	P49815
Dilution Range	IHC-P=1:100-500, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of HeLa(black line).Prima...



IHC-P analysis of rat brain using Phosph...