

Product Datasheet

Human ANGPTL7 / CDT6 Protein (orb55848)

Description	Angiopoietin-like 7 (ANGPTL7) is also known as Corneal-Derived Transcript 6 (CDT6), and is a secrete...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ANGPTL7
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.0 kDa. The protein migrates as 32-35 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_066969.1
Purity	85%
Source	Human ANGPTL7, His Tag (orb55848) is expressed from human 293 cells (HEK293). It contains AA Gln 27 - Pro 346 (Accession # O43827-1).
MW	39.0 kDa
NCBI	NP_066969.1
Expiration Date	6 months from date of receipt.