

Product Datasheet

eEF2 antibody (orb55862)

Description

eEF2 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human eEF2. Please contact us for the exact immunogen sequence. The peptide is available as orb393995.

Target

eEF2

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

53 kDa

Uniprot ID

[P13639](#)

NCBI

[006761](#), [006770](#)

Entrez


[8685](#)

Dilution Range


WB: 1-200-2000, IHC-P: 1:100-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat stom...



Immunohistochemical staining of paraffin...



Immunofluorescence image of rat stomach

...