


Product Datasheet

Guanylate Cyclase beta1 subunit antibody (orb55865)

Description

Guanylate Cyclase beta1 subunit antibody



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat


Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications



IHC-P staining of rat cerebellum tissue ...

Immunogen

KLH conjugated synthetic peptide derived from human Guanylate Cyclase beta1 subunit. Please contact us for the exact immunogen sequence. The peptide is available as orb393998.

Target


Guanylate Cyclase beta1 subunit

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml



Immunohistochemical staining of paraffin...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype


IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography



IHC-P image of rat cerebellum tissue
usi...

MW

53 kDa

Uniprot ID

[Q02153](#)

NCBI

[006761](#), [006770](#)

Entrez

[8685](#)

Dilution Range

WB: 1-200-2000, IHC-P: 1:100-500, IF/ICC: 1/100-500

Expiration Date

12 months from date of receipt.