

Product Datasheet

14-3-3 Theta (5G1) antibody (orb559067)

Description	14-3-3 Theta (5G1) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant human 14-3-3 Theta protein, full length.
Target	14-3-3 Theta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	28 kDa
Uniprot ID	P27348
Expiration Date	12 months from date of receipt.