



Product Datasheet

Cytokeratin 4 antibody (orb559075)


Description	Cytokeratin 4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Cytokeratin 4 protein (350-480aa)
Target	Cytokeratin 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	57 kDa
Uniprot ID	P19013
Dilution Range	WB=1:500-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
esophageal
using...



IHC-P
analysis of
mouse
esophageal
using...



IHC-P
analysis of
Rat tongue
using
Cytok...