

Product Datasheet

Cyclin B2 antibody (orb559088)

Description	Cyclin B2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Recombinant human Cyclin B2 protein (1-150a)
Target	Cyclin B2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	45 kDa
Uniprot ID	O95067
Expiration Date	12 months from date of receipt.