



Product Datasheet

Histone H2A.Z antibody (orb559098)

Description	Histone H2A.Z antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Histone H2A.Z protein (1-100aa)
Target	Histone H2A.Z
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	13 kDa
Uniprot ID	POC055
Dilution Range	WB=1:500-1000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human lung
carcinoma
U...



Western
blot: Lane
1: Human
HeLa
cell; Lane...