

Product Datasheet

HIF-1 Beta antibody (orb559100)

Description	HIF-1 Beta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HIF-1 beta
Target	HIF-1 beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	87 kDa
Uniprot ID	P27540
Expiration Date	12 months from date of receipt.