

Product Datasheet

Histone H2A.X (phospho-Ser139) antibody (orb559106)

Description

Histone H2A.X (phospho-Ser139) antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human H2AX around the phosphorylation site of Ser139 (QA(p-S)QE)

Target

Phospho-H2AX (Ser139)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Recombinant

MW

16 kDa

Uniprot ID

P16104

Dilution Range


WB=1:500-2000, IHC-P=1:50-200

Expiration Date


12 months from date of receipt.



IHC-P
analysis of
mouse
stomach
using Ph...



IHC-P
analysis of
rat testis
using
Phosp...



Western
blot: Lane
1: MCF-7
(Human)
Cell; L...