



Product Datasheet

C-Myc (phospho-Ser62) antibody (orb559118)

Description	C-Myc (phospho-Ser62) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human C-Myc around the phosphorylation site of Ser62 (PL(p-S)PS)
Target	phospho-C-Myc (Ser62)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	49 kDa
Uniprot ID	P01106
Dilution Range	IHC-P=1:50-200, WB=1:500-2000, IHC-F=1:50-200, ICC=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat testis
using
phospho...



IHC-P
analysis of
human
laryngeal
carcin...