



Product Datasheet

c-Jun (phospho-Ser63) antibody (orb559120)

Description	c-Jun (phospho-Ser63) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Jun around the phosphorylation site of Ser63 (LT(p-S)PD)
Target	phospho-c-Jun (Ser63)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	36 kDa
Uniprot ID	P05412
Dilution Range	WB=1:500-2000, IHC-P=1:50-100
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
esophageal
cance...



Western
blot: Lane
1: NIH/3T3
(Mouse)
Cell...