



Product Datasheet

RPS6KA1 (phospho-Ser380) antibody (orb559128)

Description	RPS6KA1 (phospho-Ser380) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPS6KA1 around the phosphorylation site of Ser380 (GF(p-S)FV)
Target	Phospho-RPS6KA1 (Ser380)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	81 kDa
Uniprot ID	Q15418
Dilution Range	WB=1:1000-5000, IHC-P=1:100-500, ICC=1:50
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat testis
using
Phosp...



IHC-P
analysis of
human
colon using
Phos...