

Product Datasheet

RPS6KA1 (phospho-Thr359 + Ser363) antibody (orb559129)

Description	RPS6KA1 (phospho-Thr359 + Ser363) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPS6KA1 around the phosphorylation site of Thr359/Ser363 (SR(p-T)PKD(p-S)PG)
Target	Phospho-RPS6KA1 (Thr359 + Ser363)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	81 kDa
Uniprot ID	Q15418
Expiration Date	12 months from date of receipt.